H250











The H250 Convertible Chair is a no-lift, multi-functional device that helpscaregivers to more efficiently and effectively performs af epatient handling practices that may help to improve patient satisfaction, safety and dignity, and reduce pain. Its Tilt-n-Space capability and therapeutic cushioning allows increased blood circulation helping to prevent skin breakdown and other complications of immobility.

Unlimited Repositioning

The H250 Convertible Chair is a new type of transfer and repositioning device that easily converts from a stretcher to an upright mobile chair in seconds. It allows a single caregiver the ability to move a patient from a bed to a chair with little effort. Its Tilt-n-Space feature allows it to move positioning on the frame thus allowing the patient to remain stationery without the angle of the hips, knees or ankles afftected. Providing a means for weight shifts and position changes in unlimited positions reduces the patient's risk of skin breakdown by redistributing pressure. Since the H250 Convertible Chair does not require a sling, there is no pressure on existing wounds or interference with healing.

Easy Transfers

The H250 Convertible chair when used with the patient transfer system allows supine transfers in four easy steps (see illustrations). Its compact overall size makes it ideal for home use and its exceptional maneuvrability makes it easy to move from room to room.

Quality

Our products are designed and manufactured to meet all the requirements of relevant standards. At every stage of the development process, every single product is given extensive attention in design, development and manufacturing. This means that a user of this product feels confident with its reliability and the caregiver feels confident with its safety.

HUMAN CARE SWEDEN (HQ)

Årstaängsvägen 21B 117 43 Stockholm Phone:- +46 8 665 35 00 Fax:- +46 8 665 35 10 info.se@humancaregroup.com

HUMAN CARE AUSTRALIA

Unit 4, 49 Bertie Street Port Melbourne, VIC 3207 Phone: +61 1300 132 736 Fax: +61 [03] 9645 4954 info.au@humancaregroup.com

HUMAN CARE CANADA

10-155 Colonnade Road Ottawa, ON K2E 7K1 Phone: +613.723.6734 Fax: +613.723.1058 info.ca@humancaregroup.com

HUMAN CARE NEDERLAND

Elspeterweg 124 8076 PA Vierhouten Phone: +31 577 412 171 Fax: +31 577 412 170 info.nl@humancaregroup.com

HUMAN CARE UNITED STATES

8006 Cameron Road, Suite K Austin, TX 78754 Phone: +512.476.7199 Fax: +512.476.7190 info.us@humancaregroup.com



PTS System

Patient transfer system used for transfer to/from a bed to/from the chair. Part no. b2004e



Transfer Sheets

Standard or bariatric sizing. Part no. b1410 (standard), b1420 (bariatric).



Castors

6" Directional or Lock Replacement Part no.: t4071 (Directional), t4070 (Lock)



Manufactured to all relevant standards, including

Transfers



Place the chair by the bed in the stretcher position.



Step 2

Attach the bed sheet to the transfer bar.



Transfer the patient by turning the handle.



Incline the chair to transport patient.

Features



Headrest

Adjustable in height. depth and angle for proper support.



Tilt-n-Space

Designed for pressure relief and positioning. Moves position without affecting body angle of patient.



Adjustable Leg/Footrest Adjustable to 90 ° angle

and 6" for added length and comfort.



Anti-Bacterial Fabric

Anti-bacterial, antimicrobial, breathable material that can be easily cleaned or disinfected in seconds.

Technical Specifications

H250			
Part No.	B1002-BL	Overall Height (Prone)	68cm/27"
Seat Width	43cm/17"	Overall Height (Wheelchair)	135cm - 152cm / 53" - 60"
Overall Width	61cm/24"	Back Angle	0° to 70°
Overall Length (Prone)	184cm - 203cm / 72.5" - 80"	Tilting Angle	0° to 17°
Overall Length (Wheelchair)	124cm - 170cm / 49" - 67"	Castors	15cm/6"
Headrest (Adjustable Height)	0cm - 17.8cm / 0" - 7"	Product Weight	62.2kg/123lb
Footrest (Adjustable Height)	0cm - 15.2cm/0" - 6"	Weight Capacity	113.4kg/250lb



Watch our chair transfer video

This product is manufactured in accordance with all relevant ISO standards.







Follow Us on Social Media













