PMF-9185

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Table of Contents

1 General Information

- 1.1 Product versions
- 1.2 Accessories
- 1.3 Scope of delivery
- 1.4 Storage, transport
- 1.5 Service life / lifecycle

2 Safety information

- 2.1 Symbols used
- 2.2 Function check

3 Using the product

- 3.1 Indication and further use
- 3.2 Contraindications
- 3.3 Application risks
- 3.4 Areas of application
- 3.5 Application illustration
- 3.6 Application
- 3.7 Approved accessories, combination with other products

4 Product marking / CE

5 Cleaning, care and disinfection

6 Technical Data

7 Warranty

8 Disposal information

1 General Information

Thank you for choosing a Petermann product.

Please read the user manual carefully before using the product for the first time!

The current understanding in modern occupational health and safety is that people in nursing professions experience significant work-related physical strain. Back injuries are common. There are many situations in which a patient has to be transferred from a bed, a wheeled stretcher, ambulance stretcher, an X-Ray table or a wheelchair to a different stretcher or seat. The strain on the user and the patient / care-dependent can be reduced significantly through the use of transfer aids.

The use of transfer aids should be conditional upon the user having learned to use the equipment properly. The user has to decide for themselves what equipment is suitable to the patient's situation and the intended sequence of movements.

The Rollboard Vision is part of a product line which is used specifically to support the repositioning and/or transfer of patients in hospitals and other health care facilities.

The Rollboard Vision is a safe transfer aid, which facilitates transfers from bed to bed or other recumbent positions in a manner pleasant and comfortable for the patient as well as simple and back-friendly for the users.

The product may only be used for the specified purpose and under the specified conditions.

The information in this user manual is for the patient who is being transferred and moved, as well as the user performing the transfer / position change.

Contact us if you require additional information or assistance using the Petermann product.

All versions of the Rollboards Vision are X-ray and MRI compatible and have successfully passed testing by TÜV SÜD.

The cover PMF-9186 has limited X-ray capability, as loops and seams may be visible in the case of multiple layers.

All Petermann Rollboards Vision are patented under the following patent no.:

- DE 10 2009 048 319
- PCT/EP2010/005877

1.1 Product versions

A

Rollboard Vision with Textile Cover:

Hygienically non-critical and semi-critical applications and locations

Simple wipe disinfection
Cover machine washable at up to 60 -80 °C

B

Rollboard Vision with PEM Disposable Cover:

Hygienically semi-critical and critical applications, application methods and locations

PEM Disposable Cover, usable without Textile Cover, for simple and fast use

Item No.:	PMF-9050 Neo Genius	Item No.:	PMF-9089
Dimensions:	850 x 450 x 10 mm, 2x foldable Set consisting of: 1 x Board PMF-9050 1 x Textile Cover PMF-9051	Dimensions:	950 x 500 x 10 mm, foldable Set consisting of: 1 x Board PMF-9085 1 x PEM Disposable Cover PMF-9090
Item No.:	PMF-9075 Neo	Item No.:	PMF-9099
Dimensions:	850 x 450 x 10 mm, foldable Set consisting of: 1 x Board PMF-9075 1 x Textile Cover PMF-9076	Dimensions:	950 x 500 x 10 mm, not foldable Set consisting of: 1 x Board PMF-9095 1 x PEM Disposable Cover PMF-9090
Item No.:	PMF-9085	Item No.:	PMF-9153
Dimensions:	950 x 500 x 10 mm, foldable Set consisting of: 1 x Board PMF-9085 1 x Textile Cover PMF-9086/96	Dimensions:	1500 x 600 x 10 mm, foldable Set consisting of: 1 x Board PMF-9150 1 x PEM Disposable Cover PMF-9152
Item No.:	PMF-9095	Item No.:	PMF-9157
Dimensions:	950 x 500 x 10 mm, not foldable Set consisting of: 1 x Board PMF-9095 1 x Textile Cover PMF-9086/96	Dimensions:	1500 x 600 x 10 mm, not foldable Set consisting of: 1 x Board PMF-9155 1 x PEM Disposable Cover PMF-9152
Item No.:	PMF-9150	Item No.:	PMF-9175
Dimensions:	1500 x 600 x 10 mm, foldable Set consisting of: 1 x Board PMF-9150 1 x Textile Cover PMF-9151	Dimensions:	1800 x 500 x 10 mm, foldable Set consisting of: 1 x Board PMF-9170 1 x PEM Disposable Cover PMF-9172

Item No.:	PMF-9155	Item No.:	PMF-9189
Dimensions:	1500 x 600 x 10 mm, not foldable	Dimensions:	1950 x 500 x 10 mm, foldable
	Set consisting of: 1 x Board PMF-9155 1 x Textile Cover PMF-9151		Set consisting of: 1 x Board PMF-9180 1 x Disposable Cover PMF-9172
Item No.:	PMF-9170		
Dimensions:	1800 x 500 x 10 mm, foldable		
	Set consisting of: 1 x Board PMF-9170 1 x Textile Cover PMF-9171		
Item No.:	PMF-9180		
Dimensions:	1950 x 500 x 10 mm, foldable		
	Set consisting of: 1 x Board PMF-9180 1 x Textile Cover PMF-9181		
Item No.:	PMF-9190		
Dimensions:	1950 x 500 x 10 mm, not foldable		
	Set consisting of: 1 x Board PMF-9190 1 x Textile Cover PMF-9181		
Item No.:	PMF-9185 with pulling loops and handles		
Dimensions Board:	1800 x 500 x 10 mm, foldable		
Dimensions Replacement cover:	1800 x 500 x 8 mm		
additional reinforcement:	1800 x 360 x 7 mm (grey)		
	Set consisting of 1 x Board PMF-9170 1 x Textile Cover PMF-9186		
	Replacement Cover PMF-9186 The cover has limited X-ray compatibility, as loops and seams may be visible in case of multiple layers!		

1.2 Accessories

C

Accessories for the Rollboard Vision with Textile

Cover:

Reordering OnlyU® PEM disposable cover:

OnlyU® Disposable Protection Foil

Item No.: PMF-0191

Dimensions per piece: 900 x 545 mm
110 pieces per roll

Suitable for: PMF-9085, PMF-9095

OnlyU® PEM Disposable Cover

Item No.: PMF-9077

Dimensions per piece: 475 x 850 mm
50 pieces folded in
dispenser pouch

Suitable for: PMF-9075

OnlyU® Disposable Protection Foil

Item No.: PMF-0192

Dimensions per piece: 1800 x 545 mm
55 pieces per roll

Suitable for: PMF-9170

OnlyU® PEM Disposable Cover

Item No.: PMF-9090

Dimensions per piece: 900 x 530 mm
120 pieces per roll

Suitable for: PMF-9085, PMF-9095

OnlyU® PEM Disposable Cover

Item No.: PMF-9152

Dimensions per piece: 1500 x 610 mm
70 pieces per roll

Suitable for: PMF-9150

OnlyU® PEM Disposable Cover

Item No.: PMF-9172

Dimensions per piece: 1800 x 530 mm
60 pieces per roll

Suitable for: PMF-9170

1.3 Scope of delivery

Option A

Rollboard Vision incl. Textile Cover

User manual

Option B

Rollboard Vision incl. PEM Disposable Cover

User manual

1.4 Storage, transport

Storage/transport		
Board	5 °C – 40 °C	protect against sunlight
Textile Cover	5 °C – 40 °C	protect against sunlight
PEM Disposable Cover	10 °C – 25 °C at 50 – 75 % relative humidity	protect against sunlight
Disposable Protection Foil	5 °C – 40 °C	protect against sunlight

1.5 Service life / lifecycle

If the product is used in accordance with the safety instructions, then the expected product service life is **five years**. However, the effective service life may vary depending on the frequency and intensity of use and the number of cleaning and disinfection cycles.

Perform regular function checks on the product and be sure to observe the safety information (see item 2). If applicable, also note the separate information for the disposable protection foils, the PEM disposable covers and/or the PEM Disposable Sliding foil with regard to the shelf life.

2 Safety information

WARNING indicates a potential danger which could lead to personal injury or material damage.

- It is the responsibility of the organisation / the user to determine, on the basis of the circumstances (weight / health condition of the patient) and the capabilities of the user, whether several users may be required for the transfer and/or whether additional safety measures are required.
- Defective products may no longer be used.
- The PEM Disposable Cover and the Protection Foil are disposable products and may only be used for one person.
- The PEM Disposable Cover and the Protection Foil may not be reprocessed. Reprocessing of the foil can lead to material damage and thus result in harm to the patient.
- The PEM Disposable Cover and the Protection Foil may not remain underneath the patient after the transfer.
- Any use other than the intended use of the PEM Disposable Cover is prohibited!
- Keep the PEM Disposable Cover and the Protection Foils away from children! Risk of suffocation.
- Risk of slipping! Ensure that the patient doesn't slide off the bed, stretcher, etc. during repositioning or transferring.
- Risk of accident! When using the Rollboard Vision, always ensure that the bed, the stretcher or any sort of mobile support whatsoever is standing securely, i.e. that brakes on it are engaged.
- Ensure that both ends of the Rollboards Vision are securely supported on the transfer surfaces. The clear distance bridged with the Rollboard Vision must not exceed 15 cm.
- It is essential to ensure that the number of users is sufficient to ensure a safe transfer in the given situation.
- Never leave the patient lying on Rollboard Vision unattended.
- The smooth surfaces significantly reduce resistance due to friction during repositioning. Check the transfer surfaces carefully to ensure they are fixed in place securely (e.g. wheelchair, bed).
- Ensure that the patient is lying securely on the Rollboard Vision and cannot slide off, especially while lifting the Rollboard Vision. (200 kg max.)
- Ensure that the transfers are performed by the users in a back-friendly manner.
- Weight limits must be observed (200 kg for tensile and brief lifting loads), optionally 300 kg with roll load.

- Don't use any pointy objects in combination with the Rollboard Vision, since this could result in damage to the product.
- Transfers under naked skin must be aided with a protective layer (e.g. towel).
- Transfer products are subject to a natural ageing process. Inspect the equipment regularly to ensure proper functionality.
- The use of non-standard cleaning agents or bleaching agents, alcohols, textile softeners or other chemical laundry additives may lead to permanent deterioration of the Rollboard Vison properties.

2.1 Symbols used

Rollboard Vision/Board

	The clear distance bridged must not exceed 15 cm		max. load capacity see item 6
	Warning		Manufacturer
	Observe the user manual		Batch no. / lot no.
	Conformity mark		Identification code
	Item number		Medical device label

Cover

	60 – 80 °C laundry		Do not tumble
	Do not iron		Do not bleach

OnlyU® Disposable Cover/Disposable Protection Foil

	For single use only		max. load capacity see item 6
	Warning		Manufacturer
	Storage temperature		Keep away from sunlight
	Moisture limit		Minimum durability
	Observe the user manual		Batch no. / lot no.
	Conformity mark		Identification code
	Item number		Medical device label

2.2 Function check

Repositioning products are also subject to a natural ageing process. Inspect the product regularly and carefully. Ensure that the product is in sound condition and that the material and the seams are not damaged. Check the product functionality by applying pressure and/or a load to the handles and ensure that the material withstands the load. Check the stability by bending the board carefully and ensuring that it is elastic, yet firm and sturdy. The smooth surface of the special sliding material allows for easy sliding and moving. Check these properties – ideally in comparison to a new product.

➤ **Defective products may no longer be used!**

3 Using the product

3.1 Indication and further use

Due to the wide range of possible applications, the indications cannot be tied to a specific clinical picture, since the Rollboards Vision are suitable for all applications involving people who require assistance during sideways transfers.

The Rollboard Vision is particularly well suited to gentle, pain-free and low-friction sideways transfers of people with restricted mobility as well as obese people. In particular, bumps or gaps which occur when two support surfaces are placed next to one another can be bridged.

The transfer of recumbent persons can be supported in emergency situations, in order to minimise patient movement when a diagnosis isn't available yet.

This product Rollboard Vision with textile Cover may be reused.

Rollboards are intended for reuse. It is not a single use product. The associated Disposable Protection Foil is a single use product and must not be reprocessed. Reprocessing the Disposable Protection Foil can lead to material damage and thus to damage to health of the patient

Before reusing the product, be sure to perform a thorough function check (see item 2.2) as well as cleaning or disinfection of the product. (see item 5)

3.2 Contraindications

Contraindications cannot be tied to a specific clinical picture. Observe the safety information under item 2.


3.3 Application risks

Untrained user

The Rollboard Vision may be lifted and/or pushed incorrectly if the user hasn't been trained in or is unfamiliar with correct use of the Rollboard Vision. In addition, movements may be performed too quickly or jerkily. This could damage the health of the patient and/or user. Observe the safety information under item 2.


3.4 Areas of application

General acute situations, rescue situations, operating theatres, X-ray departments, general stations


3.5 Application illustration

Rollboard Vision

Rollboard Vision with loops and handles

3.6 Application

3.6.1 Placement

- Place the patient on their side in a position their medical condition allows.
- Place the Rollboard Vision partway under the patient from the side.
- Place the patient on their back.

3.6.2 Transfer

- Work beside the bed.
- Push the patient onto the Rollboard Vision.
- Push the patient with Rollboard Vision in the desired transfer direction.
- If a draw sheet is used, the patient can be pulled in the desired direction on the Rollboard Vision using the draw sheet.

3.6.3 Removal

- Briefly place the patient slightly on their side.
- Remove the Rollboard from under the patient.

3.6.4 Notes

- Transfers for patients up to max. weight of 300 kg (roller load max. 300 kg)
- Pulling the Rollboard Vision by the handles during the transfer is possible up to a max. load of 200 kg (tensile load max. 200 kg).
- The Rollboard Vision can be operated by four users when transferring very obese persons.
- The Rollboard Vision can be lifted briefly during a transfer to get over small bumps under load. (200 kg max.) The product must not be used to carry patients.
- The Rollboard Vision can be pulled up or down, moved or briefly lifted out of an unfavourable transfer situation. (200 kg max.)
- Transfers are also possible with small clearances.
- Transfers are also possible if the two surfaces are at different heights.
- Transfers of patients with drains or respiratory distress syndrome are possible.

3.6.5 Additional notes for Rollboard Vision with pulling loops and handles

- Larger transfer distance
- Low contact transfer thanks to pulling loops

In order to meet the different hygiene requirements in clinical settings, there are three different versions of the Rollboard Vision which are used differently:

A. Rollboard Vision with Textile Cover

Before using the Rollboard Vision, ensure that the included Textile Cover is fitted onto the Rollboard body. Then the Rollboard Vision can be used as described (see item 3.6).

B. Rollboard Vision with PEM Disposable Cover (without Textile Cover)

The Rollboard Vision can also be used with the PEM Disposable Cover. The exceptional sliding properties of these covers simplify even difficult transfers in combination with the Rollboard Vision. The PEM Disposable Cover is made from tear-resistant and latex-free material.

The PEM Disposable Cover is shipped on a perforated roll with 60, 70 or respectively 120 pieces per roll (item no.: PMF-9172, PMF-9152, PMF-9090 or PMF-9077 – depending on the type of Rollboard Vision). The PEM Disposable Cover is shipped ready to use in a dispenser box. Open the lid of the box and turn the foil roll in the box until you see the end of the foil. Remove the tape from the foil. Pull the end of the foil out of the box. Then pull the foil up out of the box to the perforated line and then tear the foil off. Pull the PEM Disposable Cover over the Rollboard Vision. Then the Rollboard Vision can be used as described (see item 3.6).

C. Rollboard Vision with Textile Cover and Disposable Protection Foil

The Rollboard Vision can easily be used with the Disposable Protection Foil. The Disposable Protection Foil is shipped on a perforated roll (900 or respectively 1800 mm, depending on the respective Rollboard Vision size). Pull the Protection Foil over the Rollboard Cover (don't use the Disposable Protection Foil without the Rollboard Textile Cover, so as to ensure the sliding properties!). Then the Rollboard Vision can be used as described (see item 3.6).

Visit our website at www.pm-med.de or watch our comprehensive product videos on our YouTube channel.

3.7 Approved accessories, combination with other products

The Rollboard Vision can be used with the PEM Disposable Cover, the Protection Foil or the PEM-Disposable Sliding Foil according to the application requirements, also refer to the product versions under item 1.1.

Petermann transfer aids are part of a comprehensive aid product range, which are effective and back-friendly for the user and comfortable and supportive for the patients. Additional Petermann products may be used according to the situation and clinical picture. Ask our experts. We are happy to help you!

4 Product marking / CE

The Petermann Rollboard Vision is a Class I Medical Device and meets the EU-Regulation 2017/745 of the European Parliament and Council dated 05.04.2017 for medical devices, also called MDR.

Petermann GmbH is certified in accordance with DIN EN ISO 13485 by TÜV Süd, Medical Devices - Quality management systems applicable to the design and development, production and distribution of aids for people of limited mobility.

5 Cleaning, care and disinfection

Option A

Hygienically non-critical and semi-critical applications and locations

Rollboard Vision with Textile Cover
Simple wipe disinfection –
cover machine washable at up to 60 – 80 °C

5.1 General

All Petermann repositioning and transfer aids can be wet cleaned. Use warm water and a mild cleaning agent. If applicable, wipe the product off with a moist cloth and then let it dry thoroughly. Don't use any scouring or bleaching agents or any sharp-edged objects for cleaning.

When necessary, we recommend wipe disinfection with the product Trionic® (DGHM/VAH approval, further information upon request). As alternatives to our proven disinfectant Trionic®, we can also recommend the following disinfectants for wipe disinfection:

- Mikrobac® Tissues*
- Mikrobac® forte*
- Bacillol® 30 Foam
- Bacillol® 30 Tissues

When using products marked with an *, wipe down again using a water-soaked cloth after disinfection has been completed, in order to remove residues and prevent build-up. For wipe disinfection, please be sure to always follow the instructions provided by the manufacturer.

Always follow the care instructions on the label of every product as well as the information in the respective user manual.

5.2 Board

The body itself is very tough, has an inert surface and can be cleaned with conventional household cleaners (see item 5.1.). The board can be disinfected without a problem using alcohol-based disinfectant solution as well as Trionic® (DGHM/VAH approval, further information upon request). Then ensure that it dries fully. We would be happy to send you a detailed report of the board's chemical resistance upon request. Due to the temperature sensitive structure of the material, the board is not well suited to autoclave applications.

5.3 Textile Cover

The Rollboard Vision Textile Cover can be machine washed at 60° - 80° C and dried at low temperatures (45° – 50°C). Use a mild detergent! The use of bleaching agents, alcohol or fabric softeners can permanently impair the properties of the Textile Cover. Fabric materials and similar surfaces are generally incompatible with disinfectants containing alcohol and/or aldehydes. Do not autoclave!

5.4 Protection foil, PEM Disposable Cover

The foils are disposable products and may only be used for one person. They are not intended for processing or reuse and cannot be washed or treated.

Option B

Hygienically semi-critical and critical applications, application methods and locations

Rollboard Vision with PEM Disposable Cover, usable without Textile Cover, for simple and fast use

Option C

Hygienically semi-critical and critical applications and locations

Rollboard Vision with Disposable Protection Foil on top of the cover for simple and fast use

6 Technical Data

	Material	Colour	Roll-/tensile load/brief lifting load	Rollload
Rollboard Vision	Polyethylen, latex free	Graphit	200 kg	300 kg Please speak to the responsible medical device consultant if necessary
Textile Cover	Polyamide (coated), latex free	red/optional blau (different item no.)		
PEM-Disposable Cover	Polyethylene (LDPE), latex free	green	200 kg	
Disposable Protection Foil	Polyethylene (LDPE), latex free	lucent	200 kg	

Item no.	Description	Dimensions in mm	Weight	UDI
PMF-9050	Rollboard Vision Set with Textile Cover		410 g	
PMF-9050.0	Board	850 x 450 x 10	320 g	(01)04250355393517(10) x
PMF-9051	Textile Cover	800 x 450 x 1	80 g	(01)04250355390516(10) x
PMF-9075	Rollboard Vision Set with Textile Cover		610 g	
PMF-9075.0	Board	850 x 450 x 22	530 g	(01)04250355390752(10) x
PMF-9076	Cover	800 x 450 x 1	80 g	(01)04250355390769(10) x
PMF-9085	Rollboard Vision Set with Textile Cover		990 g	
PMF-9085.0	Board	950 x 500 x 22	900 g	(01)04250355390851(10) x
PMF-9086/96	Textile Cover	900 x 500 x 1	90 g	(01)04250355390868(10) x
PMF-9089	Rollboard Vision Set with PEM Disposable Protection Foil PMF-9085.0 + PMF-9090			
PMF-9090	PEM Disposable Protection Foil	108000 x 530 x 0,05	6020 g	(01)04250355390905(10) x
PMF-9095	Rollboard Vision Set with Textile Cover		990 g	
PMF-9095.0	Board	950 x 500 x 22	900 g	(01)04250355390950(10) x
PMF-9086/96	Textile Cover	900 x 500 x 1	90 g	(01)04250355390868(10) x
PMF-9099	Rollboard Vision Set with PEM Disposable Protection Foil PMF-9095.0 + PMF-9090			
PMF-9150	Rollboard Vision Set with Textile Cover			
PMF-9150.0	Board	1500 x 600 x 22	1300 g	(01)04250355391506(10) x
PMF-9151	Textile Cover	1450 x 600 x 1	120 g	(01)04250355391513(10) x

PMF-9153	Rollboard Vision Set with PEM Disposable Protection Foil PMF-9170.0 + PMF-9152			
PMF-9152	PEM Disposable Protection Foil	105000 x 610 x 0,05	6200 g	(01)04250355391520(10) x
PMF-9155	Rollboard Vision Set with Textile Cover			
PMF-9155.0	Board	1500 x 600 x 22	1300 g	(01)04250355391551(10) x
PMF-9151	Textile Cover	1450 x 600 x 1	120 g	(01)04250355391513(10) x
PMF-9157	Rollboard Vision Set with PEM Disposable Protection Foil PMF-9155.0 + PMF-9152			
PMF-9170	Rollboard Vision Set with Textile Cover			
PMF-9170.0	Board	1800 x 500 x 22	1500 g	(01)04250355391704(10) x
PMF-9171	Textile Cover	1700 x 500 x 1	150 g	(01)04250355391711(10) x
PMF-9175	Rollboard Vision Set with PEM Disposable Protection Foil PMF-9170.0 + PMF-9172			
PMF-9172	PEM Disposable Protection Foil	108000 x 530 x 0,05	6000 g	(01)04250355391728(10) x
PMF-9180	Rollboard Vision Set with Textile Cover			
PMF-9180.0	Board	1950 x 500 x 22	1800 g	(01)04250355391803(10) x
PMF-9181	Textile Cover	1800 x 500 x 1	160 g	(01)04250355391810(10) x
PMF-9185	Rollboard Vision Set with Textile Cover PMF-9170.0 + PMF-9186			
PMF-9186	Cover	1700 x 500 x 8	1000 g	(01)04250355391865(10) x
PMF-9189	Rollboard Vision Set with PEM Disposable Protection Foil PMF-9180.0 + PMF-9172			
PMF-9077	PEM Disposable Protection Foil	42500 x 475 x 0,05	1560 g	(01)04250355390776(10) x
PMF-9190	Rollboard Vision Set with Textile Cover			
PMF-9190.0	Board	1950 x 500 x 22	1800 g	(01)04250355391902(10) x
PMF-0191	Disposable Protection Foil	99000 x 545 x 0,03	3700 g	(01)04250355301918(10) x
PMF-0192	Disposable Protection Foil	99000 x 545 x 0,03	3700 g	(01)04250355301925(10) x

The manufacturing of our products is subject to highest quality standards. Due to usual manufacturing variations published dimensions and further information concerning the condition of the product are defined as target production specifications. Deviations which might occur are monitored by professional quality management procedures and do not represent a reduction in performance of the respective product.

7 Warranty

- The company Petermann warrants that the product is free from defects in material and workmanship.
- The warranty period is 3 years. **It does not cover defects attributable to abrasion, normal wear, improper use or improper cleaning.**
- The company Petermann is not liable for accidental or consequential injury or material damage.

8 Disposal information

There are no special disposal instructions for the transfer aid Rollboard Vision with Textile Cover, it can be disposed of in accordance with national and local regulations. The PEM Disposable Cover is free from hazardous substances and can be disposed of like all polyethylene materials in accordance with national and local guidelines as well as the facility standards.

Petermann GmbH
Waldweg 3
D-91601 Dombühl
Phone: +49 (0) 9868-9339-0
Email: petermann@pm-med.de

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