



The Rebel is designed for users that find it difficult to maneuver a regular rollator or require extra support. The forearm supports are adjustable in height and angle find the optimal positioning for comfort and stability.

Outstanding Stability

The Rebel RA is the ideal rollator for pain sensitive users suffering from arthritis, Parkinson's or many other neurological disorders. It features separate forearm supports and large ball-bearing mounted wheels that ensure walking is easy, even on uneven or rough surfaces. The Rebel RA can be used both indoors and outdoors and features a wide assortment of accessories including a seat and basket¹.

Unique Handles

The Rebel RA features unique front-facing, angle adjustable handles that allow users to find the optimal handle position for comfort and stability. Adjustment is easy with

integrated measurements right in the handle tube. The ergonomic design of the grips ensures even distribution of weight. The Rebel RA can also be equipped with a one-hand brake to allow a user to control hand functions on either the left or right hand side.

Adjustable Forearm Supports

The Rebel RA features very unique forearm supports that are adjustable in both height and direction. This ensures that the forearm supports can be positioned in the best possible way for user comfort and stability.

¹Standard feature in Canada

HUMAN CARE SWEDEN (HQ)

Årstaängsvägen 21B
117 43 Stockholm
Phone:- +46 8 665 35 00
Fax:- +46 8 665 35 10
info.se@humancaregroup.com

HUMAN CARE AUSTRALIA

Unit 4, 49 Bertie Street
Port Melbourne, VIC 3207
Phone : +61 1300 132 736
Fax : + 61 (03) 9645 4954
info.au@humancaregroup.com

HUMAN CARE CANADA

10-155 Colonnade Road
Ottawa, ON K2E 7K1
Phone: +613.723.6734
Fax: +613.723.1058
info.ca@humancaregroup.com

HUMAN CARE NEDERLAND

Elspeterweg 124
8076 PA Vierhouten
Phone: +31 577 412 171
Fax: +31 577 412 170
info.nl@humancaregroup.com

HUMAN CARE UNITED STATES


8006 Cameron Road, Suite K
Austin, TX 78754
Phone: +512.476.7199
Fax: +512.476.7190
info.us@humancaregroup.com

Accessories/Features



Tray

Practical for resting lightweight items. Maximum load 2kg/4.4lb. Part no. 10601



Basket

Maximum load 4.5kg/10lb. Ideal for shopping. Part no. 10621



Reflector Kit

Includes 30 self-adhesive reflectors. Part no. 91090



Slow-Down Brake

Prevents the rollator from moving too fast and enables a controlled walking speed for the user. Part no.: 10620




Light kit

Illuminate your way while walking. Part no. 91091



Security Lock

Flex cord that protects rollator from theft while unattended. Part no. 91092



Cane holder

Simplifies transport of a cane or umbrella. Part no. 10603




One-Hand Brake

Allows user to control hand functions with one hand (install left/right). Special order item. Part no. Low 10653, Standard 10654




Curb Climber

Helps raise rollator onto a curb. Part no. 10612



Oxygen Holder

Holds oxygen tube firmly in place. Part no. 10615



Seat Pad

State-of-the-art technology keeps pad cool in summer and warm in winter. Part no. 91093

Note: Not all accessories are available for sale in all countries.

This product is manufactured in accordance with all relevant ISO standards.



Technical Specifications

Description	Part No	Width between Handles	Overall Width	Handle Height (Min-Max)	Product Weight
Rebel RA 67 Low	22311	40.6cm/16"	67cm/26.4"	86cm-104cm/ 33.8" - 40.9"	11.3kg/24.9lb
	22311-1 (with seat)				
Rebel RA 67 Nor	22321	40.6cm/16"	67cm/26.4"	104cm - 122cm/ 40.9" - 48"	11.5kg/25.3lb
	22321-1 (with seat)				
Rebel RA 72 Low	22411	50.8cm/20"	72cm/28.3"	86cm-104cm/ 33.8" - 40.9"	11.4kg/25.1lb
	22411-1 (with seat)				
Rebel RA 72 Nor	22421	50.8cm/20"	72cm/28.3"	104cm - 122cm/ 40.9" - 48"	11.6kg/25.6lb
	22421-1 (with seat)				

Description (All Models)	
Seat Height	Low 54cm/21.3" Nor 62cm/24.4"
Wheel Size	20cm x 4cm/7.9" x 1.6"
Maximum User Weight	125kg / 275.6lb

Follow Us on Social Media



facebook.com/humancareworld



instagram.com/humancaregroup



linkedin.com/company/human-care-ab



pinterest.com/humancaregroup



youtube.com/channel/UC7fyInh3zTjRx0P1IDf2Da0