



The Carl-Oskar Trolley features a wide, comfortable forearm support platform table that allows pain sensitive users to comfortably rest their forearms while the front-facing brake handles assist to maintain better posture.




Outstanding Stability

The Carl-Oskar TW features a reactive frame that remains stable over virtually any surface. With extra large ball-bearing mounted wheels and a comfortable and firm table, the Carl-Oskar TW is the ideal option for pain sensitive users that find it difficult to use regular trolley walkers.

The Carl-Oskar TW can also be equipped with a one-hand brake to allow a user to control hand functions on either the left or right hand side.

Optional Accessories

The Carl-Oskar TW offers a vast assortment of accessories including oxygen tank holders that allows the walker to function for a number of users. The Carl-Oskar TW easily folds once the table is removed.



Unique Handles

The Carl-Oskar TW features unique front-facing, angle adjustable handles that allow users to find the optimal handle position for comfort and stability. Adjustment is easy with integrated measurements right in the handle tube. The ergonomic design of the grips ensures even distribution of weight.

HUMAN CARE SWEDEN (HQ)

Årstaängsvägen 21B
117 43 Stockholm
Phone:- +46 8 665 35 00
Fax:- +46 8 665 35 10
info.se@humancaregroup.com

HUMAN CARE AUSTRALIA

Unit 4, 49 Bertie Street
Port Melbourne, VIC 3207
Phone : +61 1300 132 736
Fax : + 61 (03) 9645 4954
info.au@humancaregroup.com

HUMAN CARE CANADA

10-155 Colonnade Road
Ottawa, ON K2E 7K1
Phone : +613.723.6734
Fax : +613.723.1058
info.ca@humancaregroup.com

HUMAN CARE NEDERLAND

Elspeterweg 124
8076 PA Vierhouten
Phone : +31 577 412 171
Fax: +31 577 412 170
info.nl@humancaregroup.com

HUMAN CARE UNITED STATES

8006 Cameron Road, Suite K
Austin, TX 78754
Phone: +512.476.7199
Fax: +512.476.7190
info.us@humancaregroup.com

Accessories/Features



Tray
Practical for resting lightweight items. Maximum load 2kg/4.4lb. Part no. 10601



Splash Guard
Protects legs and clothes. Part no. 10608



Reflector Kit
Includes 30 self-adhesive reflectors. Part no. 91090




Slow-Down Brake
Prevents the rollator from moving too fast and enables a controlled walking speed for the user. Part no.: 10620



Light kit
Illuminate your way while walking. Part no. 91091



Security Lock
Flex cord that protects rollator from theft while unattended. Part no. 91092




Cane holder
Simplifies transport of a cane or umbrella. Part no. 10603




One-Hand Brake
Allows user to control hand functions with one hand (installed left or right). Special order item. Part no. 10653



Curb Climber
Helps raise rollator onto a curb. Part no. 10612



Oxygen Holder
Holds oxygen tube firmly in place. Part no. 10615



Basket
Ideal for shopping. Part no. 10621

Technical Specifications

	Carl-Oskar TW Low 67	Carl-Oskar TW Low 72	Carl-Oskar TW Nor 67	Carl-Oskar TW Nor 72
Part no.	13315	13415	13325	13425
Handle Height (Min - Max)	95cm - 113cm/ 37.4" - 44.5"	95cm - 113cm/ 37.4" - 44.5"	107cm - 126cm/ 42.1" - 49.6"	107cm - 126cm/ 42.1" - 49.6"
Overall Width	67cm/26.4"	72cm/28.3"	67cm/26.4"	72cm/28.3"
Width between Handles	47cm/18.5"	52cm/20.5"	47cm/18.5"	52cm/20.5"
Distance between Rear Wheels	50.8cm/20"	52.1cm/20.5"	50.8cm/20"	52.1cm/20.5"
Wheel Size	30cm x 5cm/ 11.8" x 1.8"	30cm x 5cm/ 11.8" x 1.8"	30cm x 5cm/ 11.8" x 1.8"	30cm x 5cm/ 11.8" x 1.8"
Product Weight	14.5kg/31.9lb	14.6kg/32.2lb	14.7kg/32.4lb	14.8kg/32.6lb
Maximum User Weight	125kg/275.6lb	125kg/275.6lb	125kg/275.6lb	125kg/275.6lb

Note: Not all accessories are available for sale in all countries.

This product is manufactured in accordance with all relevant ISO standards.



Follow Us on Social Media

