



The Carl-Oskar is the perfect rollator for the active user. Its durable frame is easy to maneuver and features ergonomic handles for a comfortable grip. With extra large wheels, it is the perfect rollator for those outdoor activities, even during snowy winters.

#### **Outstanding Stability**

The Carl-Oskar is one of Human Care's best known and best-selling rollators worldwide. It is safe and reliable and features extra-large 30cm/11" wheels that allow the rollator to roll very easily yet remain stable on virtually any surface. The Carl-Oskar is the ideal option for rural areas, even on snowy surfaces.

#### **Unique Handles**

The Carl-Oskar features unique angle adjustable handles that allow users to find the optimal handle position for comfort and stability. The ergonomic design ensures even distribution of

weight for a stable grip. The Carl-Oskar can also be equipped with a one-hand brake to allow a user to control hand functions on either the left or right hand side.

#### **Reinforced Option**

The Carl-Oskar is available in different overall widths to accommodate a wide variety of users. The reinforced models for bariatric users feature a weight capacity of 275kg/606lb. All models come with a full assortment of accessories to meet most needs and can easily be folded for storage or transport.

#### **HUMAN CARE SWEDEN (HQ)**

Årstaängsvägen 21B  
117 43 Stockholm  
Phone:- +46 8 665 35 00  
Fax:- +46 8 665 35 10  
info.se@humancaregroup.com

#### **HUMAN CARE CANADA**

10-155 Colonnade Road  
Ottawa, ON K2E 7K1  
Phone: +613.723.6734  
Fax: +613.723.1058  
info.ca@humancaregroup.com

#### **HUMAN CARE NEDERLAND**

Elspeterweg 124  
8076 PA Vierhouten  
Phone: +31 577 412 171  
Fax: +31 577 412 170  
info.nl@humancaregroup.com

#### **HUMAN CARE UNITED STATES**

8006 Cameron Road, Suite K  
Austin, TX 78754  
Phone: +512.476.7199  
Fax: +512.476.7190  
info.us@humancaregroup.com

## Accessories/Features



**Tray**  
Practical for resting lightweight items. Maximum load 2kg/4.4lb.  
Part no. 10602



**Splash Guard**  
Protects legs and clothes.  
Part no. 10608



**Reflector Kit**  
Includes 30 self-adhesive reflectors.  
Part no. 91090



**Slow-Down Brake**  
Prevents the rollator from moving too fast and enables a controlled walking speed for the user. Part no.: 10624



**Light kit**  
Illuminate your way while walking. Part no. 91091



**Security Lock**  
Flex cord that protects rollator from theft while unattended.  
Part no. 91092




**Cane holder**  
Simplifies transport of a cane or umbrella.  
Part no. 10603



**One-Hand Brake**  
Allows user to control hand functions with one hand (installed left or right). Special order item. Part no. 10652



**Bag Kit**  
Includes one small (5kg capacity) and one large bag (10kg capacity). Part no. 91095 (for 62cm), 91096



**Oxygen Holder**  
Holds oxygen tube firmly in place. Part no. 10615



**Seat Pad**  
State-of-the-art technology keeps pad cool in summer and warm in winter.  
Part no. 91093



**Backstrap and Basket**  
Included with rollator in Canada. Backrest Part no. 10621, 10655 (for 72cm) Basket Part no. 10621

	Carl-Oskar 67	Carl-Oskar 72	Carl-Oskar Reinforced 67	Carl-Oskar Reinforced 72
Part no.	1135001BR	1145001BR	61325BR	61425BR
Handle Height (Min - Max)	77cm - 95cm/ 30.3" - 37.4"	77cm - 95cm/ 30.3" - 37.4"	77cm - 95cm/ 30.3" - 37.4"	77cm - 95cm/ 30.3" - 37.4"
Seat Height	64cm/25.2"	64cm/25.2"	64cm/25.2"	64cm/25.2"
Overall Width	67cm/26.4"	72cm/28.3"	67cm/26.4"	72cm/28.3"
Overall Length	82cm/32.3"	82cm/32.3"	82cm/32.3"	82cm/32.3"
Width between Handles	47cm/18.5"	52cm/20.5"	47cm/18.5"	52cm/20.5"
Distance between Rear Wheels	52.1cm/20.5"	52.1cm/20.5"	52.1cm/20.5"	52.1cm/20.5"
Wheel Size	30cm x 5cm/ 11.8" x 1.8"	30cm x 5cm/ 11.8" x 1.8"	30cm x 5cm/ 11.8" x 1.8"	30cm x 5cm/ 11.8" x 1.8"
Product Weight	11.4kg/25.1lb	11.5kg/25.3lb	13.3kg/29.3lb	13.4kg/29.5lb
Maximum User Weight	150kg/330.7lb	150kg/330.7lb	275kg/606.3lb	275kg/606.3lb