



Compact
Reactive frame
Versatile



The Rebel is a compact rollator ideal for indoor or outdoor use. Its reactive frame is stable, yet easy to maneuver. The Rebel features ergonomic handles that are angle-adjustable allowing the user to find their optimal hand position. With an easy-to-fold system, it can easily be transported.

Perfect Indoors or Outdoors

At home or on the go, the Rebel is the perfect choice. Its reactive frame combined with stable wheels allow the rollator to move smoothly on virtually any surface. Available in an assortment of widths, the Rebel offers a versatile rollator that can be used almost anywhere. The reinforced model for those more robust users features a weight capacity of 250kg/551lb.

Easy Folding System

The Rebel folds easily with a convenient folding handle built right into the seat and stands independently even in this folded position. Its compact size and

folding system make it an ideal option for transport and travel. Shown with the basket and backstrap¹, the Rebel features a full range of accessories to meet the diverse needs of each individual user.

Ergonomic Handles

In order to offer the best ergonomics, the handles feature a stable and comfortable grip that distributes the weight put on them evenly. In addition, the Rebel's handles are angle and height adjustable, so that an optimal hand position can be found for each user.

¹ Backstrap and basket are optional in all countries except for Canada

Accessories/Features



Tray

Practical for resting lightweight items. Maximum load 2kg/4.4lb. Part no. 10605



Splash Guard

Protects legs and clothes. Part no. 10608



Reflector Kit

Includes 30 self-adhesive reflectors. Part no. 91090



Slow-Down Brake

Prevents the rollator from moving too fast and enables a controlled walking speed for the user. Part no.: 10620



Light kit

Illuminate your way while walking. Part no. 91091



Security Lock

Flex cord that protects rollator from theft while unattended. Part no. 91092



Cane holder

Simplifies transport of a cane or umbrella. Part no. 10603




One-Hand Brake

Allows user to control hand functions with one hand (installed left or right). Special order item. Part no. 10652



Bag Kit

Includes one small (5kg capacity) and one large bag (10kg capacity). Part no. 91095 (for 62cm), 91096



Oxygen Holder

Holds oxygen tube firmly in place. Part no. 10615



Seat Pad

State-of-the-art technology keeps pad cool in summer and warm in winter. Part no. 91093



Backstrap and Basket

Included with rollator in Canada. Backrest Part no. 10621, 10655 (for 72cm) Basket Part no. 10651

Technical Specifications

	Rebel 67	Rebel 72	Rebel Reinforced 67	Rebel Reinforced 72
Part no.	2234101BR	2244101BR	71321BR	71421BR
Handle Height (Min - Max)	68cm - 86cm/ 26.8" - 33.8"	68cm - 86cm/ 26.8" - 33.8"	75cm - 93cm/ 29.5" - 36.6"	75cm - 93cm/ 29.5" - 36.6"
Seat Height	54cm/21.3"	54cm/21.3"	54cm/21.3"	54cm/21.3"
Overall Width	67cm/26.4"	72cm/28.3"	67cm/26.4"	72cm/28.3"
Overall Length	75cm/29.5"	75cm/29.5"	75cm/29.5"	75cm/29.5"
Width between Handles	47cm/18.5"	52cm/20.5"	47cm/18.5"	52cm/20.5"
Distance between Rear Wheels	52.1cm/20.5"	52.1cm/20.5"	52.1cm/20.5"	52.1cm/20.5"
Wheel Size	20cm x 4cm/ 7.9" x 1.6"	20cm x 4cm/ 7.9" x 1.6"	20cm x 4cm/ 7.9" x 1.6"	20cm x 4cm/ 7.9" x 1.6"
Product Weight	9.5kg/20.9lb	9.7kg/21.4lb	11.6kg/25.6lb	11.8kg/26.0lb
Maximum User Weight	140kg/309lb	140kg/309lb	250kg/551lb	250kg/551lb