



Our rail accessories and attachments are designed to make virtually every installation need possible. All mounting accessories are designed for fast and flexible installation without compromising safety.

Ceiling, Wall and Pole Attachments

Used for installations in concrete and wood joist ceilings or in conjunction with pendulums, our ceiling attachments can be installed into both straight and curved rails without drilling. Our wall attachments can be used with all our rails and are ideal when a ceiling installation is not possible or desired. When mounting in a ceiling is not possible, or if a temporary solution is required, we offer poles to accommodate this installation. All attachments are made of a high quality, durable aluminum.

Pendulums and Distance Straps

Pendulums are suitable for acoustic or uneven ceilings. They offer an aesthetic

solution that ensures the rail looks level and neat. The number of pendulums used depends on the spans permitted and the ceiling material. Our pendulums come in two different models in assorted lengths to accommodate most every installation need. Distance straps allow for an extension of the distance between the rail and the lift. They come in sizes of up to 100cm/39.3”.

Rail Rollers

We offer a wide variety of rail rollers including rail rollers for rail systems with in-rail charging. Our rail rollers come in different widths and designs and can be used with all our lift and rail systems.

HUMAN CARE SWEDEN (HQ)

Årstaängsvägen 21B
117 43 Stockholm
Phone: +46 8 665 35 00
Fax: +46 8 665 35 10
info.se@humancaregroup.com

HUMAN CARE AUSTRALIA

Unit 4, 49 Bertie Street
Port Melbourne, VIC 3207
Phone : +61 1300 132 736
Fax : + 61 (03) 9645 4954
info.au@humancaregroup.com

HUMAN CARE CANADA

10-155 Colonnade Road
Ottawa, ON K2E 7K1
Phone: +613.723.6734
Fax: +613.723.1058
info.ca@humancaregroup.com

HUMAN CARE NEDERLAND

Elspeterweg 124
8076 PA Vierhouten
Phone: +31 577 412 171
Fax: +31 577 412 170
info.nl@humancaregroup.com

HUMAN CARE UNITED STATES

8006 Cameron Road, Suite K
Austin, TX 78754
Phone: +512.476.7199
Fax: +512.476.7190
info.us@humancaregroup.com

Accessories



Rail Rollers

Assorted models for portable or stationary lifts. Model no. 20107 shown.



Distance Straps

Extends the distance between the rail and the lift. Assorted lengths.




Altair 150, 220 or 300
Compatible with all Altair models. Part no. 55200, 55210 or 55220.



Traverse Rollers

Assorted models for use with traverse system. Model no. 50170 shown.



Wall Attachment

When ceiling attachment is not possible, wall attachments can be used with all rails. Model no. 53025 shown.




Roomer S
Compatible with Roomer S. Part no. 55300



Pendulums

Assorted models and heights to accommodate installation in acoustic or uneven ceilings. Model no. 50190 shown.




Poles

Ideal for temporary use or when modification to the building structure is not permitted. Model no. 53030 shown (pole).




HeliQ
Compatible with all HeliQ models. Part no. 55100, 55110 or 55120



Attachments

Attachment kits for mounting rails onto various materials. Model no. 50167 and 50168 shown.



In-Rail Charging

Allows continuous charging of lift while in rail. Assorted components required.



HeliQ Power Traverse
Compatible with all HeliQ PT models. Part no. 55103, 55140 or 55150

Technical Specifications

Code	Description						
Rail Rollers							
20107	With Loop, Portable Lifts	Model no. 20107	Model no. 20529	Model no. 55618	Model no. 54300	Model no. 71353	Model no. 54302
20529	With Hook, Portable Lifts						
55618	Standard, Stationary Lifts						
54300	For 1" Gap, Stationary Lifts						
71353	In-rail charging, Stationary Lifts						
54302	In-rail charging, 1" Gap, Stationary Lifts						
Rail Rollers with Different Widths							
71381	With Loop, Portable Lifts, 31mm wide		71382	With Hook, Portable Lifts, 31 mm wide			
71384	With Loop, Portable Lifts, 37mm wide		71385	With Hook, Portable Lifts, 37 mm wide			
71387	With Loop, Portable Lifts, 39mm wide		71388	With Hook, Portable Lifts, 39 mm wide			
71930	With Loop, Portable Lifts, 41mm wide		71391	With Hook, Portable Lifts, 41 mm wide			

(Continued. on next page)

This product is manufactured in accordance with all relevant ISO standards.






















Learn more about this product



Follow Us on Social Media



Technical Specifications (Rail Accessories)

Code	Description				
Rail Rollers with Different Widths					
71383	Stationary Lifts, 31mm wide		71394	Stationary Lifts with In-Rail Charging, 31 mm wide	
71386	Stationary Lifts, 37mm wide		71395	Stationary Lifts with In-Rail Charging, 37 mm wide	
71389	Stationary Lifts, 39mm wide		71396	Stationary Lifts with In-Rail Charging, 39 mm wide	
71993	Stationary Lifts, 41mm wide		71397	Stationary Lifts with In-Rail Charging, 41 mm wide	
Traverse Rollers					
50170	Standard (Pair)	Model no. 50170	Model no. 53022	Model no. 54308	Model no. 71837
53022	Raised Rail (Pair)				
54308	Standard with In-Rail Charging (Pair)				
71837	Low (Pair)				
Rail Roller Accessories					
55654	Trolley Support Kit	Model no. 55654	Model no. 71326	Model no. 71838	
71326	Charging Station (IRC Rail Rollers)				
71838	Charging Station (Power Traverse)				
Attachments					
50167	Ceiling Attachment	Model no. 50167	Model no. 53008	Model no. 53025	Model no. 50499
53008	Wood Ceiling Attachment				
53025	Wall Attachment				
50499	Joint Attachment				
Poles					
53030	Complete Pole, 2 Pack (270cm)	Model no. 53030	Model no. 20228	Model no. 20230	
20228	Support, 2 Pack (302cm)				
20230	Support, 2 Pack (470cm)				
Pendulums					
50190	Adjustable 10cm - 14cm / 3.9" - 5.5"	Model no. 5019x	50760-360	Adj. 26cm - 36cm / 10.2" - 14.2"	Model no. 50760-xxx
50191	Adjustable 14cm - 18cm / 5.5" - 7"		50760-520	Adj. 34cm - 52cm / 13.3" - 20.5"	
50192	Adjustable 18cm - 22cm / 7" - 8.7"		50760-1400	Adj. 52cm - 140cm / 20.5" - 55.1"	
50193	Adjustable 22cm - 26cm / 8.7" - 10.2"		30125	Fixed - specify length required	
50194	Adjustable 26cm - 30cm / 10.2" - 11.8"				
Distance Straps					
20124-xx	Available lengths: 10cm / 3.9", 15cm/5.9", 20cm/7.8", 25cm/9.8", 30cm/11.8", 35cm/13.8", 40cm/15.7", 45cm/17.7", 50cm/19.7", 55cm/21.7", 60cm/23.6", 65cm/25.6", 70cm/27.6", 75cm/29.6", 80cm/31.5", 85cm/33.5", 90cm/35.4", 95cm/37.4", 100cm/39.3"	Strap Width 4cm/1.6"	Material Reinforced nylon thread, metal hooks		

(Continued on next page)

This product is manufactured in accordance with all relevant ISO standards.













Learn more about this product



Follow Us on Social Media



Technical Specifications (Rail Accessories)

Code	Description						
In-Rail Charging Accessories							
54324	IRC Power Rail	Model no. 54324	Model no. 54110	Model no. 54111	Model no. 54115	Model no. 54321	Model no. 50880
54110	Contactora for Power Rail						
54111	Contactora for Power Cable						
54115	IRC Entry Guide						
54321	Cable for Charger						
50880	Wall Charger						
Rail Finishes							
20069	Adjustable Bumper	Model no. 20069	Model no. 50168	Model no. 53049	Model no. 55608, 55609, 55610		
50168	Plate with Screw						
53049	Plate with Screw						
55608	End Cover - P87						
55609	End Cover - P120						
55610	End Cover - P180						
Note: Information on required screws, spacers and nuts can be found in the Rail Installation manual (99646)							

This product is manufactured in accordance with all relevant ISO standards.



Learn more about this product



Follow Us on Social Media

