




Solid Support  
Fast Installation  
Full Porfolio



Our rail accessories and attachments are designed to make virtually every installation need possible. All mounting accessories are designed for fast and flexible installation without compromising safety.

#### Ceiling, Wall and Pole Attachments

Used for installations in concrete and wood joist ceilings or in conjunction with pendulums, our ceiling attachments can be installed into both straight and curved rails without drilling. Our wall attachments can be used with all our rails and are ideal when a ceiling installation is not possible or desired. When mounting in a ceiling is not possible, or if a temporary solution is required, we offer poles or a door stand to accommodate this installation. All attachments are made of a high quality, durable aluminum.

#### Pendulums and Distance Straps

Pendulums are suitable for acoustic or uneven ceilings. They offer an aesthetic

solution that ensures the rail looks level and neat. The number of pendulums used depends on the spans permitted and the ceiling material. Our pendulums come in two different models in assorted lengths to accommodate most every installation need. Distance straps allow for an extension of the distance between the rail and the lift. They come in sizes of up to 100cm/39.3".

#### Rail Rollers

We offer a wide variety of rail rollers including rail rollers for rail systems with in-rail charging. Our rail rollers come in different widths and designs and can be used with all our lift and rail systems.

#### HUMAN CARE SWEDEN (HQ)

Årstaängsvägen 21B  
117 43 Stockholm  
Phone:- +46 8 665 35 00  
Fax:- +46 8 665 35 10  
info.se@humancaregroup.com

#### HUMAN CARE CANADA

10-155 Colonnade Road  
Ottawa, ON K2E 7K1  
Phone: +613.723.6734  
Fax: +613.723.1058  
info.ca@humancaregroup.com

#### HUMAN CARE NEDERLAND

Elspeterweg 124  
8076 PA Vierhouten  
Phone: +31 577 412 171  
Fax: +31 577 412 170  
info.nl@humancaregroup.com

#### HUMAN CARE UNITED STATES

8006 Cameron Road, Suite K  
Austin, TX 78754  
Phone: +512.476.7199  
Fax: +512.476.7190  
info.us@humancaregroup.com

## Accessories



### Rail Rollers

Assorted models for portable or stationary lifts. Model no. 20107 shown.




### Distance Straps

Extends the distance between the rail and the lift. Assorted lengths.



### Traverse Rollers

Assorted models for use with traverse system. Model no. 50170 shown.




### Wall Attachment

When ceiling attachment is not possible, wall attachments can be used with all rails. Model no. 53025 shown.



### Pendulums

Assorted models and heights to accommodate installation in acoustic or uneven ceilings. Model no. 50190 shown.



### Poles and Door Stands

Ideal for temporary use or when modification to the building structure is not permitted. Model no. 53030 shown (pole).



### Attachments

Attachment kits for mounting rails onto various materials. Model no. 50167 and 50168 shown.



### In-Rail Charging

Allows continuous charging of lift while in rail. Assorted components required.



### Altair 150, 220 or 300

Compatible with all Altair models. Part no. 55200, 55210 or 55220.



### Roomer S

Compatible with Roomer S. Part no. 55300



### HeliQ

Compatible with all HeliQ models. Part no. 55100, 55110 or 55120



### HeliQ Power Traverse




Compatible with all HeliQ PT models. Part no. 55103, 55140 or 55150

## Technical Specifications

Code	Description						
<b>Rail Rollers</b>							
20107	With Loop, Portable Lifts	Model no. 20107	Model no. 20529	Model no. 55618	Model no. 54300	Model no. 71353	Model no. 54302
20529	With Hook, Portable Lifts						
55618	Standard, Stationary Lifts						
54300	For 1" Gap, Stationary Lifts						
71353	In-rail charging, Stationary Lifts						
54302	In-rail charging, 1" Gap, Stationary Lifts						
<b>Rail Rollers with Different Widths</b>							
71381	With Loop, Portable Lifts, 31mm wide		71382	With Hook, Portable Lifts, 31 mm wide			
71384	With Loop, Portable Lifts, 37 mm wide		71385	With Hook, Portable Lifts, 37 mm wide			
71387	With Loop, Portable Lifts, 39 mm wide		71388	With Hook, Portable Lifts, 39 mm wide			
71930	With Loop, Portable Lifts, 41 mm wide		71391	With Hook, Portable Lifts, 41 mm wide			






(Continued on next page)

# Technical Specifications (Rail Accessories)

Code	Description						
<b>Rail Rollers with Different Widths</b>							
71383	Stationary Lifts, 31mm wide		71394	Stationary Lifts with In-Rail Charging, 31 mm wide			
71386	Stationary Lifts, 37mm wide		71395	Stationary Lifts with In-Rail Charging, 37 mm wide			
71389	Stationary Lifts, 39mm wide		71396	Stationary Lifts with In-Rail Charging, 39 mm wide			
71993	Stationary Lifts, 41mm wide		71397	Stationary Lifts with In-Rail Charging, 41 mm wide			
<b>Traverse Rollers</b>							
50170	Standard (Pair)	Model no. 50170	Model no. 53022	Model no. 54308	Model no. 71837		
53022	Raised Rail (Pair)						
54308	Standard with In-Rail Charging (Pair)						
71837	Low (Pair)						
<b>Rail Roller Accessories</b>							
55654	Trolley Support Kit	Model no. 55654	Model no. 71326	Model no. 71838			
71326	Charging Station (IRC Rail Rollers)						
71838	Charging Station (Power Traverse)						
<b>Attachments</b>							
50167	Ceiling Attachment	Model no. 50167	Model no. 53008	Model no. 53025	Model no. 50499		
53008	Wood Ceiling Attachment						
53025	Wall Attachment						
50499	Joint Attachment						
<b>Poles and Door Frame</b>							
53030	Complete Pole, 2 Pack (270cm)	Model no. 53030	Model no. 20228	Model no. 20230	Model no. 53029		
20228	Support, 2 Pack (302cm)						
20230	Support, 2 Pack (470cm)						
53029	Complete Door Frame						
<b>Pendulums</b>							
50190	Adjustable 10cm - 14cm / 3.9" - 5.5"	Model no. 5019x	50760-360	Adj. 26cm - 36cm / 10.2" - 14.2"	Model no. 50760-xxx		
50191	Adjustable 14cm - 18cm / 5.5" - 7"		50760-520	Adj. 34cm - 52cm / 13.3" - 20.5"			
50192	Adjustable 18cm - 22cm / 7" - 8.7"		50760-1400	Adj. 52cm - 140cm / 20.5" - 55.1"			
50193	Adjustable 22cm - 26cm / 8.7" - 10.2"		30125	Fixed - specify length required			
50194	Adjustable 26cm - 30cm / 10.2" - 11.8"						
<b>Distance Straps</b>							
20124-xx	Available lengths: 10cm / 3.9", 15cm/5.9", 20cm/7.8", 25cm/9.8", 30cm/11.8", 35cm/13.8", 40cm/15.7", 45cm/17.7", 50cm/19.7", 55cm/21.7", 60cm/23.6", 65cm/25.6", 70cm/27.6", 75cm/29.6", 80cm/31.5", 85cm/33.5", 90cm/35.4", 95cm/37.4", 100cm/39.3"	Strap Width 4cm/1.6"	Material Reinforced nylon thread, metal hooks				

(Continued on next page)

# Technical Specifications (Rail Accessories)

Code	Description						
<b>In-Rail Charging Accessories</b>							
54324	IRC Power Rail	Model no. 54324	Model no. 54110	Model no. 54111	Model no. 54115	Model no. 54321	Model no. 50880
54110	Contact for Power Rail						
54111	Contact for Power Cable						
54115	IRC Entry Guide						
54321	Cable for Charger						
50880	Wall Charger						
<b>Rail Finishes</b>							
20069	Adjustable Bumper	Model no. 20069	Model no. 50168	Model no. 53049	Model no. 55608, 55609, 55610		
50168	Plate with Screw						
53049	Plate with Screw						
55608	End Cover - P87						
55609	End Cover - P120						
55610	End Cover - P180						
Note: Information on required screws, spacers and nuts can be found in the Rail Installation manual (99646)							