

Mobi Pro Flexi



Portable
Three capacities
Lightweight



The Mobi Pro Flexi is an ergonomic and compact mobile assistive mobility device ideal for people with low leg strength and some trunk balance who experience difficulty walking. The Mobi Pro Flexi is an excellent safe patient handling solution for early mobilization. It offers support while standing and is an efficient device for patient transfers, improving safety and efficiency.

Rotating Seat

The Mobi Pro Flexi is an efficient and safe tool for the healthcare provider to move a patient, improving safety and caregiver time management. It features rotating padded seats and knee support for easy application. The feet are placed on a solid, firm footplate and then the patient pulls himself/herself up to a standing position using the integrated solid bracket. By simply rotating the seat cushion, the Mobi Pro becomes a comfortable seat that the patient can rest and be moved on.

Versatile Design

Using the Mobi Pro Flexi not only encourages mobilization and activity but gives the caregiver fast, stable and safe transport of a

patient to the bed, chair or the toilet and back if needed. The Mobi Pro Flexi is equipped with a mechanical mechanism to spread the legs, ensuring it will work with most seating (like a toilet or wheelchair). The vast adjustment range is easily set to three settings with a lever.

Quality

Our products are designed and manufactured to meet all the requirements of relevant standards. At every stage of the development process, every single product is given extensive attention in design, development and manufacturing. This means that a user of this product feels confident with its reliability and the caregiver feels confident with its safety.



Walking Aids



Healthcare Beds



Lifting Solutions



Hygiene



Convertible chairs

Accessories



Foot Plate

Included with device.
Durable ABS plastic.
Replacement Part no.
1590035



Castors

Double wheel 7.5cm/3".
Front Part no. 1690060
Rear with Brake Part no.
1690061



Ergonomic Grips

Made of rubber. Sold in
pairs. Replacement Part
no. 1510054

Lifting Situations



Use to/from a bed

Use to raise/lower a
patient from/to a bed.




Use to/from a chair

Use to raise/lower a
patient from/to a chair.



Use to/from the toilet

Use to raise/lower a
patient from/to a toilet.




Tested to all relevant
standards, including ISO
14971, EN12182, MDD
93/42/EEG.

Features




Easy Movement

The Mobi Pro Flexi is
easily moved by a single
caregiver. A small turning
radius makes it ideal for
smaller spaces.



Adjustable Base

Manual adjustment of
legs allows it to fit under
most seating.



Padded Rotating Seats

Allows device to be
used for seating as an
alternative to a
wheelchair.



Padded Knee Pad

Soft, cushioned knee pad
to add patient comfort
and reduce device
pressure at the knees
while standing.

Technical Specifications

Mobi Pro Flexi			
Part No.	92500	Adjustable Leg (Min - Max)	43cm - 70cm / 17" - 27.5"
Width	58.5cm / 23"	Handle Height (Inclined at Angle)	92cm - 100cm / 36.2" - 39.4"
Length	97cm / 38.2"	Weight	19.6kg / 43.2lb
Height	104cm / 41"	Weight Capacity	150kg / 330.6lb
Wheel Diameter	7.5cm / 3"		



*We improve freedom of movement
and quality of care for people.*

www.humancaregroup.com



HUMAN CARE SWEDEN (HO)
Årstaängsvägen 21B
117 43 Stockholm
+46 8 665 35 00
info@humancaregroup.com

HUMAN CARE GERMANY
Industriestraße 13, Gebäude 1/9
63755 Alzenau
+49 (0) 6023 994 936 8
infode@humancaregroup.com

HUMAN CARE CANADA
10-155 Colonnade Road
ON K2E 7K1 Ottawa
+1 613 723 6734
info.ca@humancaregroup.com

HUMAN CARE THE NETHERLANDS
Elspeterweg 124
8076 PA Vierhouten
+31 577 41 21 71
info.nl@humancaregroup.com

HUMAN CARE UNITED STATES
8006 Cameron rd, Suite K
TX 78754 Austin
+1 5 124 76 7199
info.us@humancaregroup.com