



The Carl-Oskar is the perfect rollator for the active user. Its durable frame is easy to maneuver and features ergonomic handles for a comfortable grip. With extra large wheels, it is the perfect rollator for those outdoor activities, even during snowy winters.



Outstanding Stability

The Carl-Oskar is one of Human Care's best known and best-selling rollators worldwide. It is safe and reliable and features extra-large 30cm/11" wheels that allow the rollator to roll very easily yet remain stable on virtually any surface. The Carl-Oskar is the ideal option for rural areas, even on snowy surfaces.

weight for a stable grip. The Carl-Oskar can also be equipped with a one-hand brake to allow a user to control hand functions on either the left or right hand side.

Reinforced Option

The Carl-Oskar is available in different overall widths to accommodate a wide variety of users. The reinforced models for bariatric users feature a weight capacity of 275kg/606lb. All models come with a full assortment of accessories to meet most needs and can easily be folded for storage or transport.



Unique Handles


The Carl-Oskar features unique angle adjustable handles that allow users to find the optimal handle position for comfort and stability. The ergonomic design ensures even distribution of

Accessories/Features



Tray

Practical for resting lightweight items. Maximum load 2kg/4.4lb. Part no. 10602



Splash Guard

Protects legs and clothes. Part no. 10608



Reflector Kit

Includes 30 self-adhesive reflectors. Part no. 91090



Slow-Down Brake

Prevents the rollator from moving too fast and enables a controlled walking speed for the user. Part no.: 10620




Light kit

Illuminate your way while walking. Part no. 91091



Security Lock

Flex cord that protects rollator from theft while unattended. Part no. 91092



Cane holder

Simplifies transport of a cane or umbrella. Part no. 10603




One-Hand Brake

Allows user to control hand functions with one hand (installed left or right). Special order item. Part no. 10652




Bag Kit

Includes one small (5kg capacity) and one large bag (10kg capacity). Part no. 91095 (for 62cm), 91096



Oxygen Holder

Holds oxygen tube firmly in place. Part no. 10615



Seat Pad

State-of-the-art technology keeps pad cool in summer and warm in winter. Part no. 91093



Backstrap and Basket

Included with rollator in Canada. Backrest Part no. 10621, 10655 (for 72cm) Basket Part no. 10621

	Carl-Oskar 67	Carl-Oskar 72	Carl-Oskar Reinforced 67	Carl-Oskar Reinforced 72
Part no.	1135001BR	1145001BR	61325BR	61425BR
Handle Height (Min - Max)	77cm - 95cm/ 30.3" - 37.4"	77cm - 95cm/ 30.3" - 37.4"	77cm - 95cm/ 30.3" - 37.4"	77cm - 95cm/ 30.3" - 37.4"
Seat Height	64cm/25.2"	64cm/25.2"	64cm/25.2"	64cm/25.2"
Overall Width	67cm/26.4"	72cm/28.3"	67cm/26.4"	72cm/28.3"
Overall Length	82cm/32.3"	82cm/32.3"	82cm/32.3"	82cm/32.3"
Width between Handles	47cm/18.5"	52cm/20.5"	47cm/18.5"	52cm/20.5"
Distance between Rear Wheels	52.1cm/20.5"	52.1cm/20.5"	52.1cm/20.5"	52.1cm/20.5"
Wheel Size	30cm x 5cm/ 11.8" x 1.8"	30cm x 5cm/ 11.8" x 1.8"	30cm x 5cm/ 11.8" x 1.8"	30cm x 5cm/ 11.8" x 1.8"
Product Weight	11.4kg/25.1lb	11.5kg/25.3lb	13.3kg/29.3lb	13.4kg/29.5lb
Maximum User Weight	150kg/330.7lb	150kg/330.7lb	275kg/606.3lb	275kg/606.3lb