


Alpha Bed Ladder

Efficient
Independent
Raising Aid



The Alpha Bed Ladder acts as a raising aid and supports the client with a physiological pattern of movement when raising the upper part of the body and finally when reaching an upright seating position.

The loop at the end makes it possible to attach the ladder to furniture, for instance a bed head. It has 8 handles in order to perform the raising in small efficient steps.

Technical specification

Product Name	Art. nr.	Dimensions	Material	Max lifting weight
Alpha Bed Ladder	PM-6090	170 x 8 cm 66,9" x 3,1"	Polyester	200kg / 440lb



CE IMD = 200 kg / 440lb

PETERMANN
part of the Human Care Group

These products are distributed by
Human Care AB for Petermann GmbH.



Visit www.humancaregroup.com for the latest product information and contact details to national offices and partners.