

Alpha Chair Sliding Mats

One-way glide
Easy
Comfortable



Alpha Chair Sliding Mat with fixation loops, open model / PM-2076



Alpha Chair Sliding Mat with 4 hand grips, open model / PM-2078




Alpha Chair Sliding Mat, tubular model / PM-2070




The Alpha Chair Sliding Mat is a safe transfer aid for clients who have slid forward in the (wheel)chair. The Chair Sliding Mat can be used to return patients to a stable sitting position. Depending on the circumstances, the client can place themselves in a stable position independently or with assistance.

All Chair Sliding Mats only slide in one direction and are placed between the client and the seat and remain in the chair during use. The sewn-on arrow indicates the direction of movement.



Technical specification

Product Name	Art. nr.	Dimensions	Material	Max lifting weight
Chair Sliding Mat	PM-2070	50 x 45 x cm 19,7" x 17,7"	Cotton/natural thickening, latex free	200 kg 440 lb
Chair Sliding Mat with fixation loops	PM-2076	130" x 45" cm 51,2" x 17,7"	Cotton/natural thickening, latex free	200 kg 440 lb
Chair Sliding Mat with 4 hand grips	PM-2078	135 x 43 cm 53,1" x 16,9"	Cotton/natural thickening, latex free	200 kg 440 lb



PETERMANN
part of the Human Care Group

These products are distributed by
Human Care AB for Petermann GmbH.



Visit www.humancaregroup.com for the latest product information and contact details to national offices and partners.