


Alpha X-Ray Cassette Placer

Alpha X-Ray Protection Foil

For easy placement of X-Ray cassettes



The Alpha X-Ray Cassette Placers are the perfect aid for safe and gentle handling of X-ray cassettes under the client.

The innovative principle of the rotating X-Ray Cassette Placer, where the X-Ray Cassette Placer is simply slid under the client, is a significant improvement in the use of X-ray cassettes on clients.

In order to meet the strict hygiene requirements, the Petermann X-Ray Cassette Placer can easily be disinfected. We also offer disposable Protection Foil for the X-ray cassette to facilitate a safe and efficient workflow.



Technical specification Alpha X-Ray Cassette Placer

Product Name	Art. nr.	Dimensions	Material	Max lifting weight
Alpha X-Ray Cassette Placer	PMF-8260	41 x 53,5 cm 16,1" x 21,1"	Polyamide, latex free	200 kg/ 440 lb
Alpha X-Ray Cassette Placer XL	PMF-8360	47 x 53,5 cm 18,5" x 21,1"	Polyamide, latex free	200 kg / 440 lb

    = 200 kg / 440lb

Technical specification Alpha X-Ray Protection Foil

Product Name	Art. nr.	Dimensions	Material
Alpha X-Ray Protection Foil	PMF-8261	45 x 80 cm / 17,7" x 31,4"	LDPE
Alpha X-Ray Protection Foil XL	PMF-8361	51 x 80 cm / 20" x 31,5"	LDPE

100 piece/roll



These products are distributed by Human Care AB for Petermann GmbH.

Visit www.humancaregroup.com for the latest product information and contact details to national offices and partners.