



ADP Approved  
ADP # AA0001895.

Easy to adjust  
and fold  
Modern design



## A smart and attractive walking companion

The Crossover is the latest contribution to our range of rollators. It is equipped with intelligent features and comes in a stylish design developed for individuals looking for an active lifestyle, supported by a stable rollator.



The Crossover folds very easily, making it ideal when traveling or if you're going to store it in a narrow space. Just pull the red strap on the back of the seat, and the folding is in process. The Crossover soft mesh basket also folds easily together with the rollator. If the basket is full of items, it can quickly be removed before folding. It is just as easy to remove it as to attach it.

means that cars will easily spot you when crossing the street. If you want to sit down for a rest, the wide padded backrest gives comfortable support. The seat can also be equipped with extra padding if you use the Crossover for sitting during more extended periods.



It also comes with a smart curb climber placed on either the left or right side. With this in place, you can lift the rollator's front with minimal effort, making it easy to get it up on a curb. As some extra safety precaution, we have also added some reflectors close to the front wheels, which

By pressing the red buttons on the frame, you can easily adjust the Crossover's height to set it to your individual needs.

We hope you will enjoy this exceptional walking aid as much as we do. With every detail in mind, we are sure it will be a stable companion to you during your walks.



## Features



### Padded Backrest

When sitting on the seat, the padded backrest gives comfortable support for the back.



### Curb Climber

The curb climber makes it easy to get the rollator over obstacles such as curbs.



### Anatomic Brake Handles

The anatomic brake handles are easy and comfortable to use.



### Soft mesh basket

Foldable basket, perfect for storing items when out and walking.

## Accessories



### Foam Cushioned Pads

The pads are placed on the seat to give a more comfortable and softer surface to sit on. Part. no 90232



### Cane Holder

Simplifies transport of a cane or umbrella. The cane holder is very easy to attach. Part. no. 90230



### Slow-Down Brake

Prevents the rollator from moving too fast and enables a controlled walking speed. Part. no 90231



### Identification Tag Green

Attaches to the rollator in seconds. Part no. 90393.



### Identification Tag Canada

Attaches to the rollator in seconds. Part no. 90392.

## Technical specifications

Part Name	Part no.	Handle Height (Min - Max)	Seat Height	Overall Width	Width between Handles	Wheel Size	Product Weight
Crossover 21	90222	86-98.5 cm 33.8" - 38.8.5"	53 cm 21"	66 cm 26"	48 cm 19"	20 cm 7,8"	8 kg 17lb

Maximum User Weight: 136 kg / 300 lb



*We improve freedom of movement  
and quality of care for people.*



### HUMAN CARE SWEDEN (HQ)

Årstaängsvägen 21B  
117 43 Stockholm  
+46 8 665 35 00  
info@humancaregroup.com

### HUMAN CARE CANADA

10-155 Colonnade Road  
ON K2E 7K1 Ottawa  
+1 613 723 6734  
info.ca@humancaregroup.com

### HUMAN CARE THE NETHERLANDS

Els peterweg 124  
8076 PA Vierhouten  
+31 577 41 21 71  
info.nl@humancaregroup.com

### HUMAN CARE UNITED STATES

8006 Cameron rd, Suite K  
TX 78754 Austin  
+1 5 124 76 7199  
info.us@humancaregroup.com