



# Alpha Rollboard Vision

Gentle  
Low Friction  
Hygienic



The Rollboard Vision is particularly well suited to gentle, pain-free and low-friction sideways transfers of people with restricted mobility as well as obese people. In particular, bumps or gaps which occur when two support surfaces are placed next to one another can be bridged. The transfer of recumbent persons can be supported in emergency situations, in order to minimise patient movement when a diagnosis isn't available yet.

The Alpha Rollboard Vision impresses with its outstanding ease of handling and uncompromising safety during transfer and care, thanks to its unique combination of hygienic features, quality, form and function. The Alpha Rollboard Vision combines size and value in its ability to handle the average size and weight of clients. The Rollboard requires a cover for use, for instance the Textile Cover, but we do recommend adding any of our covers and foils to this product.

During the development of the Alpha Rollboard Vision the demands on products for hospitals, clinics and other health facilities were taken into account along with the professional experience of people in caring professions.

All of the Rollboards have been tested and approved for use in X-ray and MRI scans.


  
**PETERMANN**  
part of the Human Care Group

These products are distributed by  
Human Care AB for Petermann GmbH.



# Various sizes of Alpha Rollboard Vision

----- where the rollboard folds



## Technical specification

### Rollboard (without cover)

	Measurement	Article number
Alpha Rollboard Vision XS, double foldable	85 x 45 cm / 33,5 x 17,7"	PMF-9050.0
Alpha Rollboard Vision S, foldable	85 x 45 cm / 33,5 x 17,7"	PMF-9075.0
Alpha Rollboard Vision M, foldable	95 x 50 cm / 37,4 x 19,7"	PMF-9085.0
Alpha Rollboard Vision M, not foldable	95 x 50 cm / 37,4 x 19,7"	PMF-9095.0
Alpha Rollboard Vision foldable (Wide)	150 x 60 cm / 59 x 23,6"	PMF-9150.0
Alpha Rollboard Vision (Wide)	150 x 60 cm / 59 x 23,6"	PMF-9155.0
Alpha Rollboard Vision L, foldable	180 x 50 cm / 70,9 x 19,7"	PMF-9170.0
Alpha Rollboard Vision XL, foldable	195 x 50 cm / 76,7 x 19,7"	PMF-9180.0

### Rollboard ProCover: Textile

	Measurement	Article number
Alpha Rollboard Vision Textile ProCover M	90 x 50 cm / 35,4 x 19,7"	PMF-9084
Alpha Rollboard Vision Textile ProCover Wide	145 x 60 cm / 57,1 x 23,6"	PMF-9154
Alpha Rollboard Vision Textile ProCover L	170 x 50 cm / 66,9 x 19,7"	PMF-9174
Alpha Rollboard Vision Textile ProCover XL	180 x 50 cm / 70,9 x 19,7"	PMF-9184

### Rollboard Cover: Textile

Alpha Rollboard Vision Textile Cover XS	80 x 45 cm / 31,5 x 17,7"	PMF-9051
Alpha Rollboard Vision Textile Cover S	80 x 45 cm / 31,5 x 17,7"	PMF-9076
Alpha Rollboard Vision Textile Cover M	90 x 50 cm / 35,4 x 19,7"	PMF-9086/96
Alpha Rollboard Vision Textile Cover (Wide)	145 x 60 cm / 57,1 x 23,6"	PMF-9151
Alpha Rollboard Vision Textile Cover L	170 x 50 cm / 66,9 x 19,7"	PMF-9171
Alpha Rollboard Vision Textile Cover XL	180 x 50 cm / 70,9 x 19,7"	PMF-9181
Alpha Rollboard Cover, with the integrated pull loops, L	170 x 50 cm / 66,9 x 19,7"	PMF-9186

### OnlyU Disposable Cover (alternative cover)

OnlyU Disposable Cover S (50 pieces in dispenser)	47,5 x 85 cm / 18,7 x 33,5"	PMF-9077
OnlyU Disposable Cover M (120 pieces per roll)	90 x 53 cm / 35,4 x 20,9"	PMF-9090
OnlyU Disposable Cover (Wide) (70 pieces per roll)	160 x 61 cm / 63 x 24"	PMF-9152
OnlyU Disposable Cover L & XL (60 pieces per roll)	180 x 53 cm / 70,9 x 20,9"	PMF-9172

### OnlyU Protection Foil (alternative cover)

OnlyU Protection Foil XS, S (110 pieces per roll)	90 x 54,5 cm / 35,4 x 21,6"	PMF-0191
OnlyU Protection Foil M, L & XL (55 pieces per roll)	180 x 54,5 cm / 70,9 x 21,6"	PMF-0192

Rollboard:



Covers:



### OnlyU Protection Foil

used on top of textile cover.

Art nr: PMF-0191 / 90x54,5 cm  
for Rollboard Vision M

Art nr: PMF-0192 / 180x54,5 cm  
for Rollboard Vision L, XL



### OnlyU Disposable Cover

used after removing the textile cover

Art nr: PMF-9090 / 90x53 cm  
for Rollboard Vision M

Art nr ex: PMF-9172 / 180x53 cm  
for Rollboard Vision L

